

# Eneo Totale Duo Kit

## Treatment Instructions



### **Eneo Totale Blu:**

Before using the Eneo Totale "Red light" device, use the Blue light "Totale Blu" device on areas that you want to be treated. Treat for 5 minutes. After 5 minutes you may start using the Total device as instructed. The BLU helps to disinfect any dirt/ oils/ toxins stuck under the skin and will clear up acne. Most importantly, it will give the red light a clear path to work at its absolute best.



### **Eneo Totale:**

Use the device on desired areas 3 times for the first month for 5 minutes per area. Then, transition to maintaining and using the device a few times a month. Make sure to use it with any choice of serum, oil or face/body cream.