

Avologi Lumnen Laser Device



Instructions



Avologi Lumnen Laser Device

1. Begin by charging the device for at least six hours before using.
2. Once fully charged, you will use the device three times a week for the first month. After the first month transition to using the device at least once per week.
3. A kindly reminder that this device is 100 times stronger than LED Light Therapy. The device was intended to be used on the entire body. Five minutes per treatment area is enough.
4. Key areas to target: Deep Wrinkles, Scars, Sun Damage, Thread Veins, loose and sagging skin.
5. LED light therapy targets the upper layers of the skin, while non-invasive laser technology operates through bio-stimulation. This not only reaches the deepest layers of the skin but also the underlying fat and muscle tissue. This triggers a series of natural events at the base layer of the dermis, resulting in increased facial volume, hyaluronic acid synthase, and stimulation of tissue growth.
6. Overall enjoy the benefits that this Nobel Prize, Award Winning, Patent Approved Medical Device has to offer